

PŘÍLOHA Č. 1: PROGRAM NA MĚŘENÍ KALIBRAČNÍ KŘIVKY

'Vypuštění do odpadu

ExternalPump Command (?) J0R
ExternalPump Flowrate (microliter/sec) 300
ExternalPump Valve Out
ExternalPump Empty
ExternalPump Delay Until Done

Position = 3

ExternalPump Command (?) O1R

'Aspirace vody

Multiposition Valve port 1
ExternalPump Flowrate (microliter/sec) 20
ExternalPump Aspirate (microliter) 200
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace standardu

Multiposition Valve port Position
ExternalPump Flowrate (microliter/sec) 15
ExternalPump Aspirate (microliter) 300
ExternalPump Delay Until Done
Delay (sec) 5

'Aspirace vody

Multiposition Valve port 1
ExternalPump Flowrate (microliter/sec) 10
ExternalPump Aspirate (microliter) 200
ExternalPump Delay Until Done
Delay (sec) 5

'Aspirace isopropanolu pro slabé mytí

Multiposition Valve port 4
ExternalPump Flowrate (microliter/sec) 10
ExternalPump Aspirate (microliter) 500
ExternalPump Delay Until Done
Delay (sec) 5

'Aspirace vzduchu

Multiposition Valve port 6
ExternalPump Flowrate (microliter/sec) 100
ExternalPump Aspirate (microliter) 300
ExternalPump Delay Until Done
Delay (sec) 1

'Vypuštění do odpadu

ExternalPump Command (?) O9R
ExternalPump Flowrate (microliter/sec) 250
ExternalPump Empty
ExternalPump Delay Until Done

Insert File

C:\Users\Obsluha\Desktop\LipenskáZ\ProcedureCleaningSyringeFromPort6.fia

Insert

FileC:\Users\Obsluha\Desktop\LipenskáZ\ProcedureCleaningDetectorFromPort6.fia

Insert File

C:\Users\Obsluha\Desktop\LipenskáZ\ProcedureCleaningDetectorFromPort8.fia

ExternalPump Aspirate (microliter) 300

ExternalPump Delay Until Done

'Aspirace elučního činidla

ExternalPump Command (?) O1R
Multiposition Valve port 5
ExternalPump Flowrate (microliter/sec) 25
ExternalPump Aspirate (microliter) 500
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace vzduchu

Multiposition Valve port 6
ExternalPump Flowrate (microliter/sec) 100
ExternalPump Aspirate (microliter) 300
ExternalPump Delay Until Done
Delay (sec) 3

'Aspirace kyseliny

ExternalPump Command (?) J7R
ExternalPump Command (?) O2R
ExternalPump Flowrate (microliter/sec) 150
ExternalPump Aspirate (microliter) 50
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace oxidačního činidla

ExternalPump Command (?) O3R
ExternalPump Aspirate (microliter) 300
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace vzduchu

ExternalPump Command (?) O8R

ExternalPump Aspirate (microliter) 300
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace isopropanolu

ExternalPump Command (?) O6R
ExternalPump Aspirate (microliter) 200
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace indikátoru

ExternalPump Flowrate (microliter/sec) 50
ExternalPump Command (?) O4R
ExternalPump Aspirate (microliter) 100
ExternalPump Delay Until Done
Delay (sec) 1

'Aspirace pufru

ExternalPump Flowrate (microliter/sec) 100
ExternalPump Command (?) O5R
ExternalPump Aspirate (microliter) 500
ExternalPump Delay Until Done
Delay (sec) 1

'Oxidace LCV

ExternalPump Command (?) O7R
Delay (sec) 120
ExternalPump Command (?) J0R

'Měření

ExternalPump Flowrate (microliter/sec) 200
Spectrometer Analyte New Sample
Spectrometer Absorbance Scanning
ExternalPump Empty
ExternalPump Delay Until Done
Spectrometer Stop Scanning